

Regulator of G-protein signaling 1

Catalog No: #AP87176



Package Size: #AP87176-1 50ug #AP87176-2 100ug #AP87176-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Regulator of G-protein signaling 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	RGS1,B-cell activation protein BL34,Early response protein 1R20,1R20,BL34,IER1
Accession No.	Uniprot ID:Q08116
Uniprot	Q08116
GeneID	5996;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RGS1

Note: This product is for in vitro research use only