

Ribonuclease inhibitor

Catalog No: #AP87178



Package Size: #AP87178-1 50ug #AP87178-2 100ug #AP87178-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ribonuclease inhibitor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Placental ribonuclease inhibitor,Placental RNase inhibitor,Ribonuclease/angiogenin inhibitor 1,RAI,PRI,RNH
Accession No.	Uniprot ID:P13489
Uniprot	P13489
GeneID	6050;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RNH1

Note: This product is for in vitro research use only