

Prostaglandin E synthase

Catalog No: #AP87186



Package Size: #AP87186-1 50ug #AP87186-2 100ug #AP87186-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostaglandin E synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mPGES-1, Microsomal prostaglandin E synthase 1, Pges
Accession No.	Uniprot ID:Q9JM51
Uniprot	Q9JM51
GeneID	64292;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pges

Note: This product is for in vitro research use only