

## Peptidyl-prolyl cis-trans isomerase A

Catalog No: #AP87189



Package Size: #AP87189-1 50ug #AP87189-2 100ug #AP87189-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Peptidyl-prolyl cis-trans isomerase A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	PPIase A,Cyclophilin A,Cyclosporin A-binding protein,Rotamase A
Accession No.	Uniprot ID:P62936
Uniprot	P62936
GeneID	397637;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PPIA

Note: This product is for in vitro research use only