

## Peroxiredoxin-6

Catalog No: #AP87190



Package Size: #AP87190-1 50ug #AP87190-2 100ug #AP87190-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Peroxiredoxin-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	1-Cys peroxiredoxin,1-Cys PRX,Acidic calcium-independent phospholipase A2,aiPLA2,Antioxidant protein 2,Ciliary body glutathione peroxidase,Non-selenium glutathione peroxidase,NSGPx,PHGPx,AOP2,GPX,PHGPX
Accession No.	Uniprot ID:O77834
Uniprot	O77834
GeneID	282438;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PRDX6

Note: This product is for in vitro research use only