

Receptor-interacting serine/threonine-protein kinase 2

Catalog No: #AP87194



Package Size: #AP87194-1 50ug #AP87194-2 100ug #AP87194-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Receptor-interacting serine/threonine-protein kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tyrosine-protein kinase RIPK2
Accession No.	Uniprot ID:P58801
Uniprot	P58801
GeneID	192656;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ripk2

Note: This product is for in vitro research use only