

Proline dehydrogenase 1, mitochondrial

Catalog No: #AP87197



Package Size: #AP87197-1 50ug #AP87197-2 100ug #AP87197-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proline dehydrogenase 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Proline oxidase,Proline oxidase 2,p53-induced gene 6 protein,PIG6,POX2
Accession No.	Uniprot ID:O43272
Uniprot	O43272
GeneID	5625;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRODH

Note: This product is for in vitro research use only