

Pleiotrophin

Catalog No: #AP87206



Package Size: #AP87206-1 50ug #AP87206-2 100ug #AP87206-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pleiotrophin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	PTN,Pleiotrophic factor beta,PTF-beta
Accession No.	Uniprot ID:P79281
Uniprot	P79281
GeneID	397609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTN

Note: This product is for in vitro research use only