

# Melanoma antigen preferentially expressed in tumors

Catalog No: #AP87209



Package Size: #AP87209-1 50ug #AP87209-2 100ug #AP87209-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Melanoma antigen preferentially expressed in tumors
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Opa-interacting protein 4,OIP-4,Preferentially expressed antigen of melanoma,MAPE,OIP4
Accession No.	Uniprot ID:P78395
Uniprot	P78395
GeneID	23532;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PRAME

Note: This product is for in vitro research use only