

Bone marrow proteoglycan

Catalog No: #AP87211



Package Size: #AP87211-1 50ug #AP87211-2 100ug #AP87211-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bone marrow proteoglycan
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BMPG,Proteoglycan 2,Mbp-1
Accession No.	Uniprot ID:Q61878
Uniprot	Q61878
GeneID	19074;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prg2

Note: This product is for in vitro research use only