

Non-secretory ribonuclease

Catalog No: #AP87212



Package Size: #AP87212-1 50ug #AP87212-2 100ug #AP87212-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Non-secretory ribonuclease
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Eosinophil-derived neurotoxin,RNase Upl-2,Ribonuclease 2,RNase 2,Ribonuclease US,EDN,RNS2
Accession No.	Uniprot ID:P10153
Uniprot	P10153
GeneID	6036;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RNASE2

Note: This product is for in vitro research use only