

Phospholipase D2

Catalog No: #AP87214



Package Size: #AP87214-1 50ug #AP87214-2 100ug #AP87214-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase D2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PLD 2,Choline phosphatase 2,PLD1C,Phosphatidylcholine-hydrolyzing phospholipase D2
Accession No.	Uniprot ID:P97813
Uniprot	P97813
GeneID	18806;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pld2

Note: This product is for in vitro research use only