

Phospholipase D3

Catalog No: #AP87216



Package Size: #AP87216-1 50ug #AP87216-2 100ug #AP87216-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase D3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PLD 3,Choline phosphatase 3,Phosphatidylcholine-hydrolyzing phospholipase D3
Accession No.	Uniprot ID:Q2KJJ8
Uniprot	Q2KJJ8
GeneID	613932;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLD3

Note: This product is for in vitro research use only