

## Melanocyte protein PMEL

Catalog No: #AP87218



Package Size: #AP87218-1 50ug #AP87218-2 100ug #AP87218-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Melanocyte protein PMEL
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ME20-M,ME20M,Melanocyte protein Pmel 17,Melanocytes lineage-specific antigen GP100,Melanoma-associated ME20 antigen,P1,P100,Premelanosome protein,Silver locus protein homolog,D12S53E,PMEL17,SILV
Accession No.	Uniprot ID:P40967
Uniprot	P40967
GeneID	6490;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PMEL

Note: This product is for in vitro research use only