

Receptor-interacting serine/threonine-protein kinase 3

Catalog No: #AP87220



Package Size: #AP87220-1 50ug #AP87220-2 100ug #AP87220-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor-interacting serine/threonine-protein kinase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	RIP-like protein kinase 3, Receptor-interacting protein 3, RIP-3, RIP3
Accession No.	Uniprot ID:Q9Y572
Uniprot	Q9Y572
GeneID	11035;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RIPK3

Note: This product is for in vitro research use only