

Rho-associated protein kinase 1

Catalog No: #AP87228



Package Size: #AP87228-1 50ug #AP87228-2 100ug #AP87228-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Rho-associated protein kinase 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Renal carcinoma antigen NY-REN-35,Rho-associated, coiled-coil-containing protein kinase 1,Rho-associated, coiled-coil-containing protein kinase I,ROCK-I,p160 ROCK-1,p160ROCK |
| Accession No. | Uniprot ID:Q13464 |
| Uniprot | Q13464 |
| GeneID | 6093; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ROCK1

Note: This product is for in vitro research use only