

PC4 and SFRS1-interacting protein

Catalog No: #AP87233



Package Size: #AP87233-1 50ug #AP87233-2 100ug #AP87233-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PC4 and SFRS1-interacting protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lens epithelium-derived growth factor,mLEDGF,Ledgf
Accession No.	Uniprot ID:Q99JF8
Uniprot	Q99JF8
GeneID	101739;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Psip1

Note: This product is for in vitro research use only