

V(D)J recombination-activating protein 2

Catalog No: #AP87238

Package Size: #AP87238-1 50ug #AP87238-2 100ug #AP87238-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	V(D)J recombination-activating protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RAG-2,Rag-2
Accession No.	Uniprot ID:P21784
Uniprot	P21784
GeneID	19374;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rag2

Note: This product is for in vitro research use only