

# 1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase beta-1

Catalog No: #AP87240

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP87240-1 50ug #AP87240-2 100ug #AP87240-3 1mg

## Description

Product Name	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PLC-154,Phosphoinositide phospholipase C-beta-1,Phospholipase C-beta-1,PLC-beta-1,Plcb
Accession No.	Uniprot ID:Q9Z1B3
Uniprot	Q9Z1B3
GeneID	18795;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Plcb1

Note: This product is for in vitro research use only