

Radical S-adenosyl methionine domain-containing protein 2

Catalog No: #AP87244

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87244-1 50ug #AP87244-2 100ug #AP87244-3 1mg

Description

Product Name	Radical S-adenosyl methionine domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Inflammatory response gene 6 protein,Viperin,Virus inhibitory protein, endoplasmic reticulum-associated, interferon-inducible,CIG6,IRG6
Accession No.	Uniprot ID:Q9MZU4
Uniprot	Q9MZU4
GeneID	396752;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RSAD2

Note: This product is for in vitro research use only