

Retinoic acid receptor alpha

Catalog No: #AP87246



Package Size: #AP87246-1 50ug #AP87246-2 100ug #AP87246-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	RAR-alpha,Nuclear receptor subfamily 1 group B member 1,NR1B1
Accession No.	Uniprot ID:Q5FBR4
Uniprot	Q5FBR4
GeneID	480526;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RARA

Note: This product is for in vitro research use only