

GTPase HRas

Catalog No: #AP87247



Package Size: #AP87247-1 50ug #AP87247-2 100ug #AP87247-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GTPase HRas
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	H-Ras-1,Ha-Ras,Transforming protein p21,c-H-ras,p21ras,HRAS1
Accession No.	Uniprot ID:P01112
Uniprot	P01112
GeneID	3265;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HRAS

Note: This product is for in vitro research use only