

Relaxin-3

Catalog No: #AP87248



Package Size: #AP87248-1 50ug #AP87248-2 100ug #AP87248-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Relaxin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Insulin-like peptide INSL7, Insulin-like peptide 7, INSL7
Accession No.	Uniprot ID:Q8HY17
Uniprot	Q8HY17
GeneID	503836;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RLN3

Note: This product is for in vitro research use only