

Renin receptor

Catalog No: #AP87249



Package Size: #AP87249-1 50ug #AP87249-2 100ug #AP87249-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Renin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ATPase H(+)-transporting lysosomal accessory protein 2,ATPase H(+)-transporting lysosomal-interacting protein 2,Renin/prorenin receptor,Atp6ip2
Accession No.	Uniprot ID:Q9CYN9
Uniprot	Q9CYN9
GeneID	70495;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Atp6ap2

Note: This product is for in vitro research use only