

Polymerase I and transcript release factor

Catalog No: #AP87251



Package Size: #AP87251-1 50ug #AP87251-2 100ug #AP87251-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Polymerase I and transcript release factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cavin-1,cav-p60
Accession No.	Uniprot ID:P85125
Uniprot	P85125
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ptrf

Note: This product is for in vitro research use only