

# Pregnancy-specific beta-1-glycoprotein 1

Catalog No: #AP87258

Package Size: #AP87258-1 50ug #AP87258-2 100ug #AP87258-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pregnancy-specific beta-1-glycoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PS-beta-G-1,CD66 antigen-like family member F,Fetal liver non-specific cross-reactive antigen 1/2,FL-NCA-1/2,PSG95,Pregnancy-specific beta-1 glycoprotein C/D,PS-beta-C/D,CD66f,B1G1,PSBG1,PSGGA
Accession No.	Uniprot ID:P11464
Uniprot	P11464
GeneID	5669;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PSG1

Note: This product is for in vitro research use only