


Description

Product Name	QSK Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	QSK Antibody detects endogenous levels of total QSK
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human QSK
Target Name	QSK
Other Names	SIK3, KIAA0999, L19, QSK, Salt-inducible kinase 3, SIK-3, SIK family kinase 3
Accession No.	Swiss-Prot#: Q9Y2K2NCBI Gene ID:23387
Uniprot	Q9Y2K2
GeneID	23387;
Calculated MW	140kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western blot analysis of QSK expression in LO2 (A), MCF7 (B), A549 (C), CT26 (D), PC12 (E) whole cell lysates.

Note: This product is for in vitro research use only