

Ras-related C3 botulinum toxin substrate 1

Catalog No: #AP87260



Package Size: #AP87260-1 50ug #AP87260-2 100ug #AP87260-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ras-related C3 botulinum toxin substrate 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	p21-Rac1
Accession No.	Uniprot ID:P63001
Uniprot	P63001
GeneID	19353;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Rac1

Note: This product is for in vitro research use only