

RAF proto-oncogene serine/threonine-protein kinase

Catalog No: #AP87261



Package Size: #AP87261-1 50ug #AP87261-2 100ug #AP87261-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RAF proto-oncogene serine/threonine-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Proto-oncogene c-RAF,cRaf,Raf-1
Accession No.	Uniprot ID:A7E3S4
Uniprot	A7E3S4
GeneID	521196;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RAF1

Note: This product is for in vitro research use only