

Oligodendrocyte transcription factor 1

Catalog No: #AP87263



Package Size: #AP87263-1 50ug #AP87263-2 100ug #AP87263-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oligodendrocyte transcription factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Oligo1,Class B basic helix-loop-helix protein 6,bHLHb6,Class E basic helix-loop-helix protein 21,bHLHe21,BHLHB6,BHLHE21
Accession No.	Uniprot ID:Q8TAK6
Uniprot	Q8TAK6
GeneID	116448;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OLIG1

Note: This product is for in vitro research use only