

Peroxisome proliferator-activated receptor gamma

Catalog No: #AP87267



Package Size: #AP87267-1 50ug #AP87267-2 100ug #AP87267-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxisome proliferator-activated receptor gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PPAR-gamma,Nuclear receptor subfamily 1 group C member 3,Nr1c3
Accession No.	Uniprot ID:O88275
Uniprot	O88275
GeneID	25664;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pparg

Note: This product is for in vitro research use only