

Ras-related protein Rab-11A

Catalog No: #AP87269



Package Size: #AP87269-1 50ug #AP87269-2 100ug #AP87269-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Ras-related protein Rab-11A |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | Rab-11,RAB11 |
| Accession No. | Uniprot ID:P62490 |
| Uniprot | P62490 |
| GeneID | 403957; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RAB11A

Note: This product is for in vitro research use only