

Ras-related protein Rab-11A

Catalog No: #AP87270



Package Size: #AP87270-1 50ug #AP87270-2 100ug #AP87270-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ras-related protein Rab-11A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Rab-11,Rab11
Accession No.	Uniprot ID:P62492
Uniprot	P62492
GeneID	53869;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rab11a

Note: This product is for in vitro research use only