

Protein phosphatase 1 regulatory subunit 14B

Catalog No: #AP87277



Package Size: #AP87277-1 50ug #AP87277-2 100ug #AP87277-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Protein phosphatase 1 regulatory subunit 14B |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Phospholipase C-beta-3 neighbouring gene protein,PLCB3N,PNG |
| Accession No. | Uniprot ID:Q96C90 |
| Uniprot | Q96C90 |
| GeneID | 26472; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PPP1R14B

Note: This product is for in vitro research use only