

## Regenerating islet-derived protein 3-gamma

Catalog No: #AP87278



Package Size: #AP87278-1 50ug #AP87278-2 100ug #AP87278-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Regenerating islet-derived protein 3-gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	REG-3-gamma,Pancreatitis-associated protein 1B,PAP-1B,Pancreatitis-associated protein IB,PAP IB,Regenerating islet-derived protein III-gamma,REG III,PAP1B,UNQ429/PRO162
Accession No.	Uniprot ID:Q6UW15
Uniprot	Q6UW15
GeneID	130120;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:REG3G

Note: This product is for in vitro research use only