

Phosphatidylinositol-3,4,5-trisphosphate 3-phosphatase and dual-specificity protein phosphatase PTEN

Catalog No: #AP87279

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87279-1 50ug #AP87279-2 100ug #AP87279-3 1mg

Description

Product Name	Phosphatidylinositol-3,4,5-trisphosphate 3-phosphatase and dual-specificity protein phosphatase PTEN
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Mutated in multiple advanced cancers 1,Phosphatase and tensin homolog,MMAC1
Accession No.	Uniprot ID:P60483
Uniprot	P60483
GeneID	403832;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTEN

Note: This product is for in vitro research use only