

Tartrate-resistant acid phosphatase type 5

Catalog No: #AP87280



Package Size: #AP87280-1 50ug #AP87280-2 100ug #AP87280-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tartrate-resistant acid phosphatase type 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	TR-AP, Tartrate-resistant acid ATPase, TrATPase, Type 5 acid phosphatase, Uteroferrin, UF
Accession No.	Uniprot ID:P09889
Uniprot	P09889
GeneID	397414;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ACP5

Note: This product is for in vitro research use only