

GTP-binding nuclear protein Ran

Catalog No: #AP87281



Package Size: #AP87281-1 50ug #AP87281-2 100ug #AP87281-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GTP-binding nuclear protein Ran
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Androgen receptor-associated protein 24,GTPase Ran,Ras-like protein TC4,Ras-related nuclear protein,ARA24,OK/SW-cl.81
Accession No.	Uniprot ID:P62826
Uniprot	P62826
GeneID	5901;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RAN

Note: This product is for in vitro research use only