

Proteolipid protein 2

Catalog No: #AP87288



Package Size: #AP87288-1 50ug #AP87288-2 100ug #AP87288-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Proteolipid protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Differentiation-dependent protein A4,Intestinal membrane A4 protein,A4
Accession No.	Uniprot ID:Q04941
Uniprot	Q04941
GeneID	5355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLP2

Note: This product is for in vitro research use only