

## Scavenger receptor class B member 1

Catalog No: #AP87289



Package Size: #AP87289-1 50ug #AP87289-2 100ug #AP87289-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Scavenger receptor class B member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SRB1,SR-BI,Srb1
Accession No.	Uniprot ID:Q61009
Uniprot	Q61009
GeneID	20778;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Scarb1

Note: This product is for in vitro research use only