

Ras-related protein Rap-1A

Catalog No: #AP87291



Package Size: #AP87291-1 50ug #AP87291-2 100ug #AP87291-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ras-related protein Rap-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	GTP-binding protein smg p21A
Accession No.	Uniprot ID:P62833
Uniprot	P62833
GeneID	282031;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RAP1A

Note: This product is for in vitro research use only