

40S ribosomal protein S14

Catalog No: #AP87293



Package Size: #AP87293-1 50ug #AP87293-2 100ug #AP87293-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	40S ribosomal protein S14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:P62264
Uniprot	P62264
GeneID	20044;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rps14

Note: This product is for in vitro research use only