

Reactive oxygen species modulator 1

Catalog No: #AP87294



Package Size: #AP87294-1 50ug #AP87294-2 100ug #AP87294-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Reactive oxygen species modulator 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ROS modulator 1,Protein MGR2 homolog
Accession No.	Uniprot ID:P60603
Uniprot	P60603
GeneID	67067;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Romo1

Note: This product is for in vitro research use only