

Pyridoxal phosphate phosphatase

Catalog No: #AP87296



Package Size: #AP87296-1 50ug #AP87296-2 100ug #AP87296-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pyridoxal phosphate phosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PLP phosphatase,Chronophin,Cin,Plp,Plpp
Accession No.	Uniprot ID:P60487
Uniprot	P60487
GeneID	57028;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdpx

Note: This product is for in vitro research use only