

Peripheral myelin protein 22

Catalog No: #AP87297



Package Size: #AP87297-1 50ug #AP87297-2 100ug #AP87297-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peripheral myelin protein 22
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PMP-22,Growth arrest-specific protein 3,GAS-3,Gas-3,Gas3,Pmp-22
Accession No.	Uniprot ID:P16646
Uniprot	P16646
GeneID	18858;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pmp22

Note: This product is for in vitro research use only