

# Receptor-binding cancer antigen expressed on SiSo cells

Catalog No: #AP87305

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP87305-1 50ug #AP87305-2 100ug #AP87305-3 1mg

## Description

Product Name	Receptor-binding cancer antigen expressed on SiSo cells
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Cancer-associated surface antigen RCAS1,Estrogen receptor-binding fragment-associated gene 9 protein,RCAS1
Accession No.	Uniprot ID:Q865S0
Uniprot	Q865S0
GeneID	403500;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:EBAG9

Note: This product is for in vitro research use only