

60S ribosomal protein L13

Catalog No: #AP87309



Package Size: #AP87309-1 50ug #AP87309-2 100ug #AP87309-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	60S ribosomal protein L13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Breast basic conserved protein 1, Large ribosomal subunit protein eL13, BBC1, OK/SW-cl.46
Accession No.	Uniprot ID:P26373
Uniprot	P26373
GeneID	6137;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RPL13

Note: This product is for in vitro research use only