

## Runt-related transcription factor 2

Catalog No: #AP87310



Package Size: #AP87310-1 50ug #AP87310-2 100ug #AP87310-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Runt-related transcription factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Core-binding factor subunit alpha-1,CBF-alpha-1,Osteoblast-specific transcription factor 2,OSF-2,Cbfa1,Osf2
Accession No.	Uniprot ID:Q9Z2J9
Uniprot	Q9Z2J9
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Runx2

Note: This product is for in vitro research use only