

# Protein S100-B

Catalog No: #AP87311



Package Size: #AP87311-1 50ug #AP87311-2 100ug #AP87311-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein S100-B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	S-100 protein beta chain,S-100 protein subunit beta,S100 calcium-binding protein B
Accession No.	Uniprot ID:P50114
Uniprot	P50114
GeneID	20203;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:S100b

Note: This product is for in vitro research use only