

Solute carrier family 12 member 1

Catalog No: #AP87312

Package Size: #AP87312-1 50ug #AP87312-2 100ug #AP87312-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Solute carrier family 12 member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BSC1,Bumetanide-sensitive sodium-(potassium)-chloride cotransporter 2,Kidney-specific Na-K-Cl symporter,Nkcc2
Accession No.	Uniprot ID:P55014
Uniprot	P55014
GeneID	20495;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc12a1

Note: This product is for in vitro research use only